

LAB MASTER LGB5000 SERIES

5.0 Purity Nickel Plated Brass Barstock Outlet Point

Lab Master LGB5000 series outlet points are designed for use in conjunction with the LGP500 series outlet point regulators for the secondary pressure gas control of non-corrosive or toxic, high purity gases up to grade 5.0 purity (99.999%), on gas distribution systems.



- **LGB5001 Series** | Wall Mounting Outlet Point, FKM Seal
- **LGB5002 Series** | Wall Mounting Outlet Point, EPDM Seal
- **LGB5011 Series** | Bench Mounting Outlet Point, FKM Seal
- **LGB5012 Series** | Bench Mounting Outlet Point, EPDM Seal

★ FEATURES

- Isolation valve
- Non-return valve
- Wall and bench mounting
- Full range of inlet and outlet fittings

👍 SPECIFICATION

Maximum Inlet Pressure	25 bar (360 PSI)
Ports	1/4" NPT Female
Temperature Range	-20° C to +60° C (-4° F to +140° F)
Weight	0.55 kg (1.21 lbs) maximum


📄 MATERIALS

Body	Nickel Plated Brass Barstock
Seal	FKM / EPDM

*Regulator sold separately **Flow performance as complete assembly with Regulator

✓ FLOW PERFORMANCE

25 bar Inlet Pressure @ 21°C



10 bar Inlet Pressure @ 21°C

